

## Single Use Bipolar / Monopolar Forcep with Pre-attached 3m Cable

This product is disposable and is supplied sterile. Prior to the expiration date, contents are sterile if the package has not been opened or damaged.



### 1. Indications

Disposable bipolar / monopolar forceps are electrosurgical devices designed to be used in soft tissue surgical procedures



### 2. Contraindications

Any use of this instrument for tasks other than for which it is indicated will usually result in a damaged or broken instrument



### 3. Setup & Use

The cables to the surgical electrodes should be positioned in such a way that contact with PATIENT or other leads is avoided. Temporarily unused ACTIVE ELECTRODES should be stored isolated from the patient.

Connect the sterile cable to the generator only while it is in the "OFF" or "STANDBY" mode. Failure to do so may result in injury or electrical shock to the patient or operating room personnel



### 4. Caution

Due to the variability of output voltages and modes from generator to generator, do not use the disposable forcep with generators having output voltages that exceed 1100Vp-p. Refer to the appropriate electrosurgical generator manual for indications and instructions on output characteristics to ensure that all safety precautions are followed



Batch code

### 5. Warning

Connect bipolar accessories to the bipolar receptacle only and monopolar accessories to the monopolar receptacle. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions



Use by date

### 6. Power Setting Guidelines

These may vary due to differences in surgical techniques, patients, electrodes and surgical set-up. Start at the lowest power setting and increase as necessary to achieve the desired clinical effect

At the lowest power setting, test the forceps by pressing the generator's activating switch. If the generator fails to activate check the forceps connection with the cable



SCHULTZ MEDICAL  
INTERNATIONAL  
WAZIR ABAD RD, SKT-PK

### 7. Cleaning and Sterilisation

This bipolar / monopolar forcep is designed and warranted for single use only